Septiembre 2023

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Bluephage EASY KITS – ISO Kits

Easy Kits

Kit Product No.	Matrices	Volume Sample	Equivalence (Method)	Type of coliphage	No. Analysis per Kit
BPI601	Raw, treated, Wastewater, Surface water, Recreational, shellfish extracts, sediments and sludge	1mL – 10mL	ISO 10705-2 (DAL method)	Somatic	70
BPI604	Clean Water e.g. Drinking Water	100mL	Based on ISO (SAL method)	Somatic	10

Biological Materials

Kit Product No.	Description	Units
BPI603	Bacterial Host strain (WG5) + Virus bacteriophage (Φ x174)	10 + 10
BPI626	Virus bacteriophage (Φ x174)	20
BPI628	Bacterial Host strain (WG5)	20

Powder media

Kit Product No.	Description
BPI637	Modified Scholten's Broth (MSB) – 500g
BPI638	Modified Scholten's Agar (MSA) – 500g (for 550 plates approximately)


Plates with MSA layer

Kit Product No.	Description	Units
MSAI	90mm plates with MSA media	20

Reagents & materials for ISO 10705-2 + 10705-3

Kit Product No.	Description	Units
MSASSIN	Semi solid Modified Scholten's Agar prepared	10
MSB50I	Modified Scholten's Broth prepared (MSB)	10
ELUI	3mL Phage eluting solution	20

- 1.** Prewarm one Mini Erlenmeyer per sample at $(36 \pm 2)^\circ\text{C}$.
- 2.** Aseptically add the *E. coli* WQ strain tablet to the Erlenmeyer. Allow 30 min for rehydration.
- 3.** Incubate at $(36 \pm 2)^\circ\text{C}$ for 120-140 min in a shaking isolator at 100-150 rpm. After the incubation time, quickly cool the culture by placing it in a refrigerator or in melting ice. Use the inoculum culture within the same working day.
- 4.** Melt the ssMSA tubes (one tube per 7 assays), slightly loosen the cap of the tubes and put them in a boiling water bath. It is important to check that the agar is totally melted before its use. There are alternative methods to melt ssMSA tubes as the use of a microwave. However, such uses are always under client's responsibility.
- 5.** Place in a water bath at $(55 \pm 1)^\circ\text{C}$ and let them reach this temperature.
- 6.** Aseptically add 0.14 mL of the calcium chloride solution, prewarmed at room temperature.
- 7.** Aseptically add 0.23 mL of the nalidixic acid solution, previously hydrated, as described above.
- 8.** Distribute 2.5 mL aliquots into red screw-capped tubes and place in a water bath at $(55 \pm 1)^\circ\text{C}$.
- 9.** To each red screw-capped tube, add 1 mL of the original sample (or diluted or concentrated sample) prewarmed at room temperature.
- 10.** Once the inoculum culture is ready, add 1 mL to each red screw-capped tube and mix carefully avoiding the formation of air bubbles. Pour the contents onto a layer of complete MSA on a 60 mm plate prewarmed at room temperature.
- 11.** Distribute evenly, cover the plates, and allow to solidify on a horizontal, cool surface.
- 12.** Incubate the plates upside-down at $(36 \pm 2)^\circ\text{C}$ for 18±2h. Do not stack more than six plates.

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1. Preparing the sample

Prewarm the water sample at (35±2)°C. It is very important that the sample reaches this temperature before starting the analysis.

2. Preparing the host strain

Prewarm the Host Strain Vial at room temperature. Remove the seal from the cap and gently rotate it clockwise to release the host strain. Shake vigorously to resuspend the host strain and let it stand for 10 minutes, upside-down at room temperature.

3. Water sample

Measure 100 mL of the water sample. The 100 mL mark is just an orientation reference. Open the BP Testing Bottle and pour in the water sample.

4. Mixing the resuspended host strain and the sample

Open the Host Strain Vial and check that all the contents in the cap have been resuspended. Pour the host strain into the BP Testing Bottle. Close the bottle.

5. Starting the test. Resuspension of Bluephage Medium.

Rotate the BP Testing Bottle cap clockwise to make the Bluephage Medium fall inside the bottle. Shake gently to completely dissolve the Bluephage Medium in the water sample.

6. Running the test

Incubate the BP Testing Bottle at (35 ± 2)°C during 1 h. Then, transfer the BP Testing Bottle to an incubator at (30 ± 2)°C during 5 h more. Total incubation time is 6 h.

7. Reading the results

Read the test results after 6h. If the yellow colour is maintained (negative), it indicates absence of somatic coliphages. If colour changes to green-blue (positive), there is presence of somatic coliphages. The colour change always occurs before 6h if there is a presence of bacteriophages, even at very low concentrations (1 PFU / 100 mL). Do not extend the incubation and reading of tests beyond the time when the negative control presents a start of color change due to enter in the natural death phase of the culture. In case of confirmation of negative results, see Spot Test at Section 8.

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1. Preparing the sample

Prewarm the water sample at 36±2°C. It is very important that the sample reaches this temperature before starting the analysis.

2. Preparing the host strain

Prewarm the Host Strain Vial at room temperature. Press the top of the cap to ensure the internal thread fits. Rotate it clockwise to break the seal and to release the host strain. Shake vigorously to resuspend the host strain and let it stand for 10 minutes upside-down at room temperature.

3. Water sample

Measure 100 mL of the water sample. The 100 mL mark is just an orientation reference. Open the BP Testing Bottle and pour in the water sample.

4. Mixing the resuspended host strain and the sample

Open the Host Strain Vial and check that all the contents in the cap have been resuspended. Pour the host strain into the BP Testing Bottle. Close the bottle.

5. Starting the test. Resuspension of Bluephage Medium and filling of ENUMERA® Testing Tray.

Rotate the BP Testing Bottle cap clockwise to make the Bluephage Medium fall inside the bottle. Shake gently to completely dissolve the Bluephage Medium in the water sample. Carefully, pour the contents of the BP Testing Bottle in the ENUMERA® Testing Tray.

6. Sealing the ENUMERA® Testing Tray

Prewarm the ENUMERA® Sealer. Put the ENUMERA® Testing Tray in the blue mould and, carefully, seal it.

7. Running the test

Place the ENUMERA® Testing Tray on a dark mat. Start the Bluephage App, follow the instructions and take the initial picture. Incubate the ENUMERA® Testing Tray at 36 ± 2°C during 1 h, placing the membrane side up to assure gas exchange. Then, transfer the ENUMERA® Testing Tray to an incubator at 30 ± 2°C during 5.5 h more. Total incubation time is 6.5 h.

8. Reading the results

Read the test results at the end of the incubation time (6.5 h). Place the ENUMERA® Testing Tray on a dark mat and start the Bluephage App. Follow the instructions and take the final picture. The result, in PFU/100 mL, will appear on your device screen.

If the yellow colour is maintained (negative) in all wells, it indicates absence of somatic coliphages. If colour changes to green-blue (positive) in any well, there is presence of somatic coliphages. The colour change always occurs before 6.5 h if there is a presence of bacteriophages, even at very low concentrations (1 PFU / 100 mL). Do not extend the incubation and reading of tests beyond the time when the negative control presents a start of colour change due to enter in the natural death phase of the culture.

Parameter	ISO	Bluephage EASY KIT	Bluephage RAPID KIT
Sensitivity	1 PFU/1mL	1 PFU/1mL 1 PFU/100mL	15 PFU/100mL
Ease to use	*	**	***
Applicability	Water, Food, Biosolids	Water, Food, Biosolids	Water, Food, Biosolids
Availability of strains	Not included	Included	Included
Preparation of material (h) (previous steps)	40-60	0	0
Preparation of inoculum culture (h)	3-4	2	0,16
Results time (h)	18-20h	18-20	6
Hours Technician	4	1	<0.5

*Not very easy **Easy ***Very easy

Depending on the Kit (volume analysed)

Eric Queralt
CEO and co-founder

Anicet Blanch
Scientific Advisor &
Co-founder

Núria Guilera
Commercial Director

Julia Martin
CDO

Miriam Pascual
CQO

Ariadna Jorba
Biotechnologist

Mireia Azuara
Press Officer

Contacto

Núria Guilera Grandes

Commercial Director

nguiler@bluephage.com

www.bluephage.com

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